

TREAT WITH EVEN MORE CONFIDENCE


Shows either the individualized risk of distant recurrence or the individualized risk of recurrence or mortality when treated with endocrine therapy alone.

Based on a continuous scale of 0-100 and reflects individual tumor biology. The higher the patient's score, the higher the risk of distant recurrence and likelihood of chemotherapy benefit.

Contains relevant study details from either RxPONDER or SWOG 8814.

Represents the likelihood that adding chemotherapy will reduce distant recurrence risk (for applicable patients).

Helps to confirm HR+, HER2- status.



Oncotype DX Breast Recurrence Score[®] Report
Micromets & Node Positive (1-3)
Postmenopausal

EXACT
SCIENCES


Last name, First name (POST_V03-2)

Date of Birth: DD-MM-YEAR Gender: Female Report Number: OR000123456-6007 Report Date: DD-MON-YEAR
Specimen Source/ID: Breast/SP-16_0123456
Ordering Physician: Dr. First-Name I. Ordering-Physician-Lst-Name

Medical Record/Patient #: DD-MM-YEAR Client: Community Medical Center
Date of Collection: 27-Jul-2023
Specimen Received: 29-Jul-2023
Additional Recipient: Dr. First-Name I. Recipient-Physician-Last-Name
Pathologist: Dr. First-Name I. Pathologist-Last-Name

Displays the patient's individualized Recurrence Score[®] result in the context of the relevant study data from either RxPONDER or SWOG 8814, if available.

RxPONDER⁵



Methods and Limitations

The Oncotype DX Breast Recurrence Score test uses RT-PCR to provide information on prognosis and the magnitude of chemotherapy benefit to guide chemotherapy treatment decisions in patients with early-stage, HR+, HER2-negative and lymph node-negative or lymph node-positive (N1) breast cancer. Decisions on treatment should also be based on independent medical judgment of the treating physician, taking into consideration all available information concerning the patient's medical condition, in accordance with your community's standard of care.

The **Recurrence Score (RS) Result** which ranges from 0-100 is calculated from the quantitative RT-PCR analysis of 21 specific genes.

Quantitative Single-Gene Scores for quality control. The Oncotype DX test uses quantitative RT-PCR to determine the RNA expression of ER, PR and HER2, using the published validated cut-offs^{2,3}. The standard deviations of single-gene results are less than 0.5 units. The RT-PCR single-gene results may differ from ER, PR, or HER2 results reported using other methods or reported by other laboratories.

The definitions of **menopausal status** based upon the RxPONDER⁴ trial are as follows:

- **Premenopausal:** Less than 6 months since last menstrual period and no prior bilateral oophorectomy and not on estrogen replacement.
- **Postmenopausal:** Prior bilateral oophorectomy or more than 12 months since last menstrual period with no prior hysterectomy.
- If these definitions did not apply, patients were categorized as premenopausal if <50 years and postmenopausal if ≥50 years.

Distant Recurrence-Free Interval (DRFI): Time from randomization to distant recurrence or death from breast cancer.

Invasive Disease-Free Survival (IDFS): Time from randomization to invasive breast cancer recurrence (local, regional, or distant), second invasive primary cancer (breast cancer or not), or death from any cause.

References:
1. Kalinsky et al. San Antonio Breast Cancer Symposium. 2021. 2. Badve et al. J Clin Oncol. 2008. 3. Baehner et al. J Clin Oncol. 2010. 4. Kalinsky et al. NEJM. 2021. 5. Data on File

Laboratory Director(s): F. Baehner, MD; H. Bailey, MD & P. Joseph, MD

This test was developed and its performance characteristics determined by Genomic Health, Inc. It has not been cleared or approved by the FDA, nor is it currently required to be. The laboratory is regulated under CLIA and qualified to perform high-complexity testing. This test is used for clinical purposes. It should not be regarded as investigational or for research.
The Oncotype DX Breast Recurrence Score Test is an in vitro diagnostic device. CE marked under Regulation (EU) 2017/46 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices. In Japan, the test has received manufacturing and marketing approval from the Ministry of Health, Labour and Welfare as the Oncotype DX Breast Recurrence Score Program and is covered by National Health Insurance.

Genomic Health, Inc.
301 Penobscot Drive, Redwood City, CA 94063
USA/Canada +1.866.ONCOTYPE (+1.866.662.6897)

<https://precisiononcology.exactsciences.com/>
CLIA Number 05D1018272

For more information, visit OncotypeDX.com
or call 866-ONCOTYPE (662-6897).